

DATA VIEWER QUICK INSTALL INSTRUCTIONS


Prior to installing the *Data Viewer*, it is recommended that you:

- Ensure that your computer meets the system requirements.
- Close all other programs and applications.
- Ensure that you have administrative rights to install software on your computer. Consult your IT staff if you do not have administrative rights.

To Install the *Data Viewer*:


More recent versions of Windows have stricter security enhancements that only allow you to install software as read-only. In order to install and use the Data Viewer on your computer, you will have to save the software in a folder where you have administrative rights.

1. Create a new folder in your personal documents folder (or some other folder where you have administrative rights to make changes.) Name the new folder "Data Viewer". This is where you will copy the Data Viewer software.
2. Insert the Data Viewer USB/flash drive into a USB port on your computer.
3. Copy the *three* files from the flash drive and paste them to the newly created **Data Viewer** folder.




4. Copying and pasting the software to the new folder will not create a shortcut icon on your desktop. You can create a shortcut or access the database directly from the Library folder by double-clicking the **CCDViewer20d.mdb** file. This will open the software for use.

5. As newer versions of Microsoft Office have stricter security features, MS Access may require the user to Enable Content prior to using the Data Viewer software. As such, you may encounter the following error message when first using the Data Viewer.



Close the error message window by clicking on the X. You will then be able to see the security warning behind the message, and click on the **Enable Content** button.



The user can then click on the **Menu** button to load an ACF-801 data file.